Polycarbonate Capillaries User Guide

Sample preparation / PCR Mixture preparation

Samples for PCR can be handled and prepared according to the existing protocol at customer site. No changes have to be made using polycarbonate capillaries instead of the original glass capillaries from Hofmann La-Roche.

No changes have to be made preparing the PCR Mix.

The polycarbonate capillaries for PCR-reaction volumes from 10 µl to 50 µl can be used without significant changes in the PCR programme.

Handling capillaries and caps

For better/easier handling of our capillaries and caps a “Transfer-Pin” for taking the capillaries and the caps out of the rack-box is provided. The Transfer-Pin is part of our “Starter-Kit”. Use gloves to protect you PCR-reaction from contamination, e.g. DNases.

For sealing the capillaries with the caps, use the transfer pin, to transfer caps from the rack-box to the polycarbonate capillary. Ensure that each capillary is closed tightly by checking it visually.

Performing PCR / PCR programme

Our polycarbonate capillaries have a greater diameter, compared to the 20 µl glass capillaries from Roche.

To guarantee effective heat transfer from the outside of the capillaries to the middle of the PCR reaction, consider the greater diameter. It might be necessary to increase the “hold-times” at Denaturing and Annealing.

For this reason we recommend to program the Cyclor with longer hold times for the denaturing and annealing steps. We suggest to program at least 7 seconds at 95°C for the denaturing step and 5 seconds for the annealing step. Too short denaturing and annealing times will lead to a complete failure of your PCR, or at least to much higher CT-times.

If you are already using a programme with both times higher than 7 seconds, no changes have to be done, normally. Sometimes it might be better to increase the existing times by 2 seconds each.

As the effectiveness of denaturing is also depending on the GC-content of your DNA, it might be even better to prolong denaturing up to 10 seconds (only in case of no PCR signal with a denaturing time of 5 seconds).

Denaturing times at 95°C						
	glass	PC-Cap..	glass	PC-Cap.	glass	PC-Cap.
Denaturing	0 sec	10 sec	5 sec	10-12 sec	> 7 sec	> 10 sec
Annealing	0 sec	10 sec	5 sec	10-12 sec	> 7 sec	> 10 sec
Extension	10 sec	10 sec	5 sec	10-12 sec	> 7 sec	> 7 s
Denaturing time at 97°C						
	glass	PC-Cap.	glass	PC-Cap.	glass	PC-Cap.
Denaturing	0 sec	5 sec	5 sec	7-10 sec	> 7 sec	> 10 sec
Annealing	0 sec	10 sec	5 sec	10 sec	> 7 sec	> 10 sec
Extension	10 sec	10 sec	5 sec	5-8 sec	> 7 sec	> 7 s

Note: Longer denaturing times than 10-11 seconds brings no better result

Sample recovery after PCR

Some customer wants to analyse the PCR reaction on a gel. For that purpose the cap can be taken away from the capillary. The recovery of the sample is possible using pipettes with a 10 µl tip. In some cases the pipette tips are too broad (Sarstedt, AHN and other brands). In this case open the capillaries and spin the solution into a clean reaction tube.

For further info please have a look to our Web-Page: www.geneon.net or follow this link: [PC-Capillaries](#)

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Handling capillaries and caps

For better/easier handling of our capillaries and caps a "Transfer-Pin" for taking the capillaries and the caps out of the rack-box is provided. The Transfer-Pin is part of our "Starter-Kit". Use gloves to protect you PCR-reaction from contamination, e.g. DNases.

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To guarantee effective heat transfer from the outside of the capillaries to the middle of the PCR reaction, consider the greater diameter. It might be necessary to increase the "hold-times" at Denaturing and Annealing.

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