



X-Gal (5-bromo-4-chloro-3-indolyl-beta-D-galactopyranoside)

Applications:

- blue/white screening
- distinguish between positive and negative clones after transformation

Description:

X-Gal is a highly pure molecular biology grade chemical for colorimetric detection of b-galactosidase activity. It is used for "blue/white screening" of bacterial colonies to differ bacterial or phage clones after transformation. In combination with suitable bacterial cloning vectors, host strains, and IPTG, X-Gal provides an easy way to distinguish between positive and negative clones after transformation.

Quality:

Purity of >98% by HPLC. Tested in the blue/white colony screening experiment.

Note:

- Preparation of a 20 mg/ml stock solution in dimethylformamide or dimethylsulfoxide is recommended.
- Dimethylformamide dissolves some plastic materials. The direct addition of dimethylformamide containing solution to plastic Petri dishes should be avoided.

Storage:

at - 20°C in the dark; shipped at ambient temperature

Ordering information:

Cat.-no	Description	Amount
406-501	X-Gal	1 g
406-502	X-Gal	5 x 1g
406-503	X-Gal	10 x 1 g

GeneON .. a good decision ..

Contact Phone +49-(0)-621- 5720 864 Fax: +49-(0)-621-5724 462

E-Mail: <mailto:info@geneon.net> WEB: <http://www.GeneOn.net> Version: 25.10.2009 AS

Unless specified otherwise, all products of GeneON are sold for research use only.

.. a good decision ..