



ddNTPs - Dideoxynucleotide triphosphates

Cat.-No. 110-020 4x (2x500µl each ddNTP)

Used in single nucleotide polymorphisms and as 3'-end chain terminators in Sanger sequencing

Description:

ddNTPs from GeneON in combination with modified Taq DNA Polymerase are used as 3'-end chain terminators in Sanger sequencing. @ pH 8,0 in water the high purity between more than 98 % (HPLC) offers a great performance in chain termination sequencing. **Tested and designed especially for Sequencing reactions!**

Application:

2',3'-Dideoxynucleoside triphosphates inhibit the chain elongation of a given primer catalyzed by the DNA polymerase (e.g. Klenow enzyme) and are therefore used for DNA sequencing according to Sanger. Sequencing is achieved by including in each reaction a dideoxy-nucleotide that acts as a chain terminator. Four reactions are set up, each containing the same template and primer but a chain terminator specific for A, C, G or T. Because only a small amount of the chain terminator is included, incorporation into the new DNA strand is a random event. Each reaction therefore generates a collection of fragments, but every DNA strand will end at the same type of base (A, C, G or T).

SNP: multiplexing by minisequencing

Multiplexed detection of alternatively spliced transcripts

Effect on reverse transcriptase of HIV-1

Molecular dynamic calculations of DNA-polymerase beta-ddCTP complex

Spectro-Data @ pH 7,5 ±0.5, Tris-HCL

ddATP: λ_{max} 259 nm, ϵ 15.1 L mmol⁻¹ cm⁻¹

ddCTP: λ_{max} 271 nm, ϵ 8.9 L mmol⁻¹ cm⁻¹

ddGTP: λ_{max} 252 nm, ϵ 14.2 L mmol⁻¹ cm⁻¹

ddTTP: λ_{max} 267 nm, ϵ 9.6 L mmol⁻¹ cm⁻¹

Molecular Weight and Formula:

ddATP: 475.18 g/mol C₁₀H₁₆N₅O₁₁P₃ (acid free)

ddCTP: 451.16 g/mol C₉H₁₆N₃O₁₂P₃ (acid free)

ddGTP: 491.18 g/mol C₁₀H₁₆N₅O₁₂P₃ (acid free)

ddTTP: 466.17 g/mol C₁₀H₁₇N₂O₁₃P₃ (acid free)

Concentration: 10 mM – 11 mM; >98 % HPLC

Note:

- spin all reagents in the vial before opening and after setting up the Termination Mix

Quality tests:

The ddNTPs are tested:

- for the absence of DNases and RNases

- for successfully sequencing reaction

Transport: on blue ice

Storage: few days at room temperature; for long term storage @ -20 °C;

Ordering details:

Cat.-no	Description	Amount
110-020	Set of ddNTP	8 x 500µl (8 x 5 µmol)

GeneON .. a good decision ..

Contact Phone +49-(0)-621- 5720 864 Fax: +49-(0)-621-5724 462

E-Mail: <mailto:info@geneon.net> WEB: <http://www.GeneOn.net> Version: 10.2009 UK update: 19-07-2023

Unless specified otherwise, all products of GeneON are sold for research use only.

.. a good decision ..