

Random Primer Cat-No: S300, 30 µg

Description:

The primer is a mixture of single-stranded random hexanucleotides with 5'- and 3'-hydroxyl ends. Random Hexamer Primers can be used for first-strand cDNA synthesis and cloning. The primer is supplied in lyophilized form.

Quantity: 30 µg

Storage: store at - 20 °C

Ordering information

Cat.-no	Description	Amount
S300	Random Primer	30 µg

.. a good decision ..

GeneON .. a good decision ..

Contact Phone +49-(0)-621- 5720 864 Fax: +49-(0)-621-5724 462

E-Mail: <mailto:info@geneon.net> **WEB:** <http://www.taq-dna.com/> Version: 25.10.2009 AS

Unless specified otherwise, all products of GeneON are sold for research use only.