



Yeast extract powder, 500 gr  
Product.-Code: CE169

**Specification**

<b>Product</b>	<b>Yeast extract</b>
Brand	GeneON
Quality	for molecular biology
Product number	CE169, CE170
Storage	below 25°C
<b>Parameter</b>	<b>Specification</b>
DNases/RNases/Proteases	not detectable
Ash	max. 14 %
Total N	11 - 12 %
pH (1 % in water, 25°C)	6.8 - 7.2
Heavy metals (as Pb)	max. 0.0001 %
Loss on drying	max. 5 %
Coagulated protein	not detectable
Sodium chloride	max. 0.3 %
Ludwigshafen, 17.01.2019	

*.. a good decision ..*