

Product	NADH disodium salt for molecular biology
Brand	GeneON-Bioscience
Product number	CE198 ff.
Synonym(s)	β-Nicotinamide adenine dinucleotide
Formula	C ₂₁ H ₂₇ N ₇ Na ₂ O ₁₄ P ₂
Molecular weight	709.41 g/mol
CAS number	606-68-8
Storage	-20°C
Retest period (months)	36
Parameter	Specification
Assay (enzym.)	min. 95 %
Appearance	yellowish powder
Water (K.F.)	max. 8.0 %
UV Absorption	
A250/A260 (pH 10)	0.82 ± 0.03
A280/A260 (pH 10)	0.23 ± 0.02
Ludwigshafen	21. Jan 2022